



IHE Integration Statement

<i>Vendor</i>	<i>Product Name</i>	<i>Version</i>	<i>Release Date</i>
RamSoft Inc	PowerServer, Gateway	4.7	November 10, 2009

This product implements all of the transactions specified in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

<i>Integration Profiles Implemented</i>	<i>Actors Implemented</i>	<i>Options Implemented</i>
Scheduled Workflow	Image Display	None
Consistent Presentation of Images	Image Display, Print Composer	None
Access to Radiology Information	Image Display, Report Reader	None
Key Image Note	Image Display	None
Simple Image and Numeric Report	Report Creator	None
Post-Processing Workflow	Image Display	None
Reporting Workflow	Image Display	None
Portable Data for Imaging	Portable Media Creator, Image Display, Report Reader, Print Composer, Portable Media Importer	None
NM Image	Image Display	None
Mammography Image	Image Display	None
Patient Information Reconciliation	Report Creator, Image Display	None
Presentation of Grouped Procedures	Image Manager	None

RamSoft

Internet address for vendor's IHE information:

<http://software.ramsoft.biz/manuals/IHE.pdf>

Links to Standards Conformance Statements for the Implementation

HL7

<http://software.ramsoft.biz/manuals/HL7%20Conformance.pdf>

DICOM

<http://software.ramsoft.biz/manuals/DICOM%20Conformance.pdf>

Links to general information on IHE

In Europe:

<http://www.rsna.org/ihe>

In Europe:

<http://www.ihe-europe.org>

In Japan:

<http://www.jira-net.or.jp/ihe-j>