

### ***IHE Integration Statement***

<i><b>Vendor</b></i>	<i><b>Product Name</b></i>	<i><b>Version</b></i>	<i><b>Release Date</b></i>
RamSoft Inc	PowerServer, Gateway	5.0	August 12, 2011

This product implements all of the transactions specified in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

<i><b>Integration Profiles Implemented</b></i>	<i><b>Actors Implemented</b></i>	<i><b>Options Implemented</b></i>
Scheduled Workflow (SWF)	Image Manager / Image Archive	None
	Image Display	None
	Evidence Creator	None
Patient Information Reconciliation (PIR)	Image Manager / Image Archive	None
Consistent Presentation of Images (CPI)	Image Display	None
	Evidence Creator	None
	Print Composer	None
	Image Manager / Image Archive	None
Access to Radiology Information (ARI)	Image Display	None
	Image Manager / Image Archive	None
	Evidence Creator	None
Key Image Note (KIN)	Image Display	None
Simple Image and Numeric Report (SINR)	Report Creator	None
Reporting Workflow (RWF)	Evidence Creator	None
	Report Creator	None

	Image Manager / Image Archive	None
Evidence Documents (ED)	Image Display (Report Creator)	None
	Portable Media Creator	None
	Portable Media Importer	None
Portable Data for Imaging (PDI)	Image Display	None
	Report Creator	None
	Image Archive / Manager	None
	Image Display	None
NM Image (NM)	Image Archive / Manager	None
	Image Display	None
Mammography Image (MAMMO)	Image Archive / Manager	None
	Image Display	None
	Print Composer	None
Import Reconciliation Workflow (IRWF)	Image Archive / Manager	None
	Importer	Unscheduled Import
Consistent Time (CT)	Time Client	None
<b>Internet address for vendor's IHE information:</b> <a href="http://www.ramsoft.com/conformance/IHE.pdf">www.ramsoft.com/conformance/IHE.pdf</a>		
<b><i>Links to Standards Conformance Statements for the Implementation</i></b>		
<b>HL7</b>	<a href="http://www.ramsoft.com/conformance/HL7_Conformance.pdf">www.ramsoft.com/conformance/HL7_Conformance.pdf</a>	
<b>DICOM</b>	<a href="http://www.ramsoft.com/conformance/DICOMConformance.pdf">www.ramsoft.com/conformance/DICOMConformance.pdf</a>	
<b><i>Links to general information on IHE</i></b>		
In North America: <a href="http://www.rsna.org/IHE">www.rsna.org/IHE</a>	In Europe: <a href="http://www.ihe-europe.org">www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a>